


NPWT Clinical Applications



Lifotronic
— Caring for Better Life —

NPWT Dressing Kit



Diversified

Sterile

Easy-to-use

Lifotronic Technology Co., Ltd.
Address: Lifotronic Tower, No. 8, Qiuzhi East Road, Longhua District, Shenzhen, 518110, P.R.China
Tel: 86-755-29016066 Fax: 86-755-29060036 Email: inter@lifotronic.com
Web: en.lifotronic.com

Follow us on   

 (SSE: 688389)

NPWT Dressing Kit

Lifotronic NPWT Dressing Kit includes the dressing foam, film and tube. It is sterilized and for single use. Not made with natural rubber latex.It is intended to be used in conjunction with Lifotronic NPWT systems to deliver negative pressure to the wound.

An open-pored, reticulated polyurethane foam for use with NPWT therapy units, designed to promote granulation tissue formation and enhance exudate removal. Dressings can be used on a wide variety of wounds and easily trimmed to fit the contours of deep or irregularly shape wounds.

Features

Dressing Foam

- Made from polyurethane foam with excellent bio-compatibility
- Bigger aperture (400-600um)
- Good drainage effect, not easy to block




Film

- Soft and flexible to fit uneven surfaces
- Waterproof and breathable for better isolation from external pollutants



Drainage tube

- Pressure up to 500mmHg to avoid blockage
- Three kinds of specifications---single tube (drainage type), double tube (pressure type), triple tube (irrigation type)



Dressing Kit

Product code	Product description	Every sterile pouch contains
For NP -100		
PU Black Foam Dressing Kit		
03-S02A000-00001	Dressing Kit Small Thick, 10-pack	1 foam dressing (100x75x33mm) 1 sheet of film 1 set of drainage single tube 1 Y-connector
03-S02A000-00002	Dressing Kit Medium Thick, 10-pack	1 foam dressing (180x125x33mm) 2 sheets of film 1 set of drainage single tube 1 Y-connector
03-S02A000-00003	Dressing Kit Large Thick, 10-pack	1 foam dressing (260x150x33mm) 2 sheets of film 1 set of drainage single tube 1 Y-connector
For NP -600		
PU Black Foam Dressing Kit		
03-S02A000-00006	Dressing Kit Small Thick, 10-pack	1 foam dressing (100x75x33mm) 1 sheet of film 1 set of drainage double tube 1 Y-connector
03-S02A000-00007	Dressing Kit Medium Thick, 10-pack	1 foam dressing (180x125x33mm) 2 sheets of film 1 set of drainage double tube 1 Y-connector
03-S02A000-00008	Dressing Kit Large Thick, 10-pack	1 foam dressing (260x150x33mm) 2 sheets of film 1 set of drainage double tube 1 Y-connector
03-S02A000-00009	Dressing Kit Medium Thin, 10-pack	1 pre-cut foam dressing (180x125x16mm) 2 sheets of film 1 set of drainage double tube 1 Y-connector
03-S02A000-00010	Dressing Kit Large Thin, 10-pack	1 pre-cut foam dressing (260x150x16mm) 2 sheets of film 1 set of drainage double tube 1 Y-connector

Dressing Kit

Product code	Product description	Every sterile pouch contains
For NP -800		
PU Black Foam Dressing Kit		
03-S02A000-00011	Dressing Kit Small Thick, 10-pack	1 foam dressing (100x75x33mm) 1 sheet of film 1 set of drainage triple tube 1 Y-connector
03-S02A000-00012	Dressing Kit Medium Thick, 10-pack	1 foam dressing (180x125x33mm) 2 sheets of film 1 set of drainage triple tube 1 Y-connector
03-S02A000-00013	Dressing Kit Large Thick, 10-pack	1 foam dressing (260x150x33mm) 2 sheets of film 1 set of drainage triple tube 1 Y-connector
03-S02A000-00014	Dressing Kit Medium Thin, 10-pack	1 pre-cut foam dressing (180x125x16mm) 2 sheets of film 1 set of drainage triple tube 1 Y-connector
03-S02A000-00015	Dressing Kit Large Thin, 10-pack	1 pre-cut foam dressing (260x150x16mm) 2 sheets of film 1 set of drainage triple tube 1 Y-connector



Dressing Foam

Product code	Product description	Every sterile pouch contains
PU Black Foam		
03-S02A000-00078	Dressing Small -Thick Foam only, 20-pack	1 foam dressing (100x75x33mm)
03-S02A000-00079	Dressing Medium -Thick Foam only, 20-pack	1 foam dressing (180x125x33mm)
03-S02A000-00081	Dressing Large -Thick Foam only, 10-pack	1 foam dressing (260x150x33mm)
03-S02A000-00080	Dressing Medium -Thin Foam only, 20-pack	1 foam dressing (180x125x16mm)
03-S02A000-00082	Dressing Large -Thin Foam only, 10-pack	1 foam dressing (260x150x16mm)
12-S02A000-00060	Dressing,Medium,Foam only	1 PU Black Foam Medium Size (180x125x16mm)

Film and Tube

Product code	Product description	Every sterile pouch contains
03-S02A000-00064	Film only, 100-pack	1 sheet of film (300*300mm)
03-S02A000-00077	Single tube only, For NP-100, 20-pack	1 set of drainage single tube
03-S02A000-00076	Double tube only, For NP-600, 20-pack	1 set of drainage double tube
03-S02A000-00075	Triple tube only, For NP-800, 20-pack	1 set of drainage triple tube
09-S02A000-00060	Film only,for Abdominal Foam	4 drapes (300x200mm)
11-S02A000-00060	Film only,for Abdominal Foam	1 organ protective film
10-S02A000-00060	Tube only,for Abdominal Foam	1 double tube drainage type